














Sof-Lex™ Finishing and Polishing System

Abrasive Grades	Sof-Lex	Sof-Lex Extra Thin (XT)	Sof-Lex Finishing Strip	Sof-Lex Polishing Brush
Superfine	 1958SF 1982SF 1981SF	 2382SF 2381SF	SF  1956	 1751
Fine	 1958F 1982F 1981F	 2382F 2381F	F  1956	
Medium	 1958M 1982M 1981M	 2382SM 2381M	M  1954	
Coarse	 1958C 1982C 1981C	 2382C 2381C	C  1954N 1954	

Sof-Lex Finishing Brush

Intro Kit 1750

- 12 Sof-Lex Finishing Brushes
- 2 Sof-Lex Finishing Brush Mandrels – Right Angle
- 30 ea. Assorted Sof-Lex Extra Thin (XT) Contouring and Polishing Discs – Coarse, Medium, Fine and Superfine – 1/2"
- 1 Sof-Lex Disc Mandrel – Right Angle



Mandrel Refill 1755 – Right Angle

- 3 Sof-Lex Finishing Brush Mandrels

Mandrel Refill 1756 – Friction Grip

- 3 Sof-Lex Finishing Brush Mandrels

Brush Refill 1751

- 24 Sof-Lex Finishing Brushes



3M ESPE

3M ESPE AG · ESPE Platz
82229 Seefeld · Germany
E-mail: info3mespe@mmm.com
Internet: <http://www.3mespe.com>

3M, ESPE and Sof-Lex are trademarks of 3M or 3M ESPE AG. Used under license in Canada.

3M ESPE

From Front to Back to
Everywhere in Between



Sof-Lex

Finishing and Polishing System


3M ESPE is a recognized leader in dental innovation. Introduced in the mid-1970s, the Sof-Lex Finishing and Polishing System was a revolutionary breakthrough – providing dentists with the ultimate advanced restoration finishing solution.

With a color-coded, sequential system of discs, strips and easy-to-use mandrels, the Sof-Lex system continues to provide unparalleled, efficiently achieved results.

Now, the renaissance continues...by pushing innovation to the tip with the Sof-Lex Finishing Brush.

Introducing the Sof-Lex™ Finishing Brush from 3M ESPE


The Sof-Lex Finishing Brush combines the once separate finishing and polishing procedures of posterior composite restorations into one, easy step. No stopping to switch shapes or grits saves time and effort – making restorative procedures much more efficient for you.



The brush bristles contain aluminum oxide abrasives for highly effective polishing. The soft, elastomer bristles are gentle to gingival tissue and conform as they travel over concave and convex surface anatomies to deliver a uniform finish comparable to Sof-Lex Discs.



The Sof-Lex Finishing Brush is also autoclavable. As part of the 3M ESPE Sof-Lex Finishing and Polishing System, the finishing brush delivers the superior results you demand – in less time and with less effort than you expect.



Sof-Lex Finishing and Polishing Discs and Strips

Sof-Lex Discs

- Reversible Discs**
 Makes aligning disc easier for polishing various surfaces.
- Fast and Easy Disc Change**
 Simply snaps securely on or off mandrel with no need for alignment.
- Color-Coded**
 Sequenced from dark (coarse) to light (superfine) for easily followed step-by-step process.
- Two Sizes**
 For a variety of needs.
 13 mm or 9 mm
- Maximized Abrasive Surface**
 Mandrel is low profile. Disc has a small round eyelet.
- Choice of disc Thickness and Flexibility**
 Applicable to a variety of uses with regular or extra thin construction.



Pick up the pace without sacrificing great results.

With coarse to superfine grits clearly color-coded for ease of procedure sequencing and quick, pop-on/off disc interchangeability, Sof-Lex finishing and polishing strips and discs deliver exceptional, time-proven results with time-saving advantages.

Sof-Lex Finishing and Polishing Discs




Sof-Lex Extra Thin Finishing and Polishing Discs



Sof-Lex Strips

For finishing and polishing interproximal areas, Sof-Lex strips round-out the system.

- Easy to Insert**
 Abrasive is center gapped on the strip for easy interproximal insertion without abrading contacts.
- Two Grits per Strip**
 Each strip performs two functions to save time.
- Color-Coded**
 Sequenced from dark to light for easily followed step-by-step process.
- Choice of Grits and Widths**
 A choice of four grits and two widths.
- Resists Tearing**
 Polyester backing resists tearing while retaining flexibility and being gentle on gingiva.



Sof-Lex Finishing and Polishing Strips

