Shaded PMMA Disc Temporary Crown and Bridge Provisorische Kronen und Brücken Couronne et bridge temporaire Temporanea per corone e ponti Temporal de corona y puente Tijdelijke kroon en brug Dočasná korunka a můstek Laikinoji karūnėlė ir tiltelis Coroană și punte provizorii Dočasná korunka a mostík

Instructions for Use
Gebrauchsanweisung
Mode d'emploi
Istruzioni per l'uso
Instrucciones de uso
Gebruiksaanwijzing
Návod k použití
Naudojimo instrukcijos
Instrucțiuni de utilizare
Návod na používanie

THE DENTAL
SOLUTIONS
COMPANY™

Intended purpose

For temporary crowns and bridges for dental prosthetic restorations.

en

Indications for Use

Dentsply Sirona Shaded PMMA Discs are indicated for manufacturing temporary crowns and bridges by a CAD/CAM system.

Contraindications

1. May not be used for patients with hypersensitivity to methacrylate or any of the other ingredients.

Intended user

Dental technicians

Intended Patient Population and medical conditions

The product is marketed to dental labs and multiple patient specific restorations can be created. This product is neither single use nor reusable.

Clinical Benefit/Patient specific outcome

Provisional restoration provides patient with the ability to masticate.

Warnings

1. Possible cross-reactions or interactions of this material with other material(s) already present in the oral environment must be taken into consideration by the dentist.

Precautions

1. Keep product dust away from eyes.
2. After use, wash your hands.
3. Do not swallow the product.
4. Do not inhale dust particles during grinding.
5. Use local vacuum suction and suitable mouth/face protection during manual machining at the workplace.

Adverse Reactions

1. If properly processed and used, adverse effects are highly unlikely. However, reactions of the immune system (such as allergies) to the substances contained in the material or localized paraesthesia (such as taste disturbances or irritation of the oral mucosa) cannot be completely ruled out.
2. Residual dust during milling operation may cause irritation.

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations. For a summary of the safety and clinical performance of this product, see <https://ec.europa.eu/tools/eudamed>

Technical Specifications

Meets applicable requirements of ISO 10477: 2020[E] Dentistry - Polymer-based Crown and Veneering Materials.

Composition

PMMA Shaded 99,6% Polymethylmethacrylate cross-linked polymer with 0,4% pigments and fluorescing agent.

Requirements for Milling Temporary Crowns & Bridges

- > Millable in any dental mill that accepts a 98.5mm disc and has milling options for PMMA materials
- > Minimum thickness: 1mm occlusion and 0.6mm cervical
- > Minimum connection diameters: Anterior: 3.6mm; Posterior: 4.4mm
- > The provisional can be fixed with any cement that has an indication in the IFU for the cementation of provisional restorations.

Step-by-Step Instructions for Milling Temporary Crowns & Bridges

A. Procedure for Fabrication

1. Place the disc on the machine following the equipment manufacturer's instructions.
2. After milling, remove the disc from the equipment and remove the sprues.
3. Finish and polish using conventional techniques.

Disposal

Dispose in accordance with local and national regulations.

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Dentsply LLC
York, PA 17404 U.S.A.
DeguDent GmbH