

CZR Liquids

Každá komponenta produktu musí být smíchána se správným liquidem, jak je uvedeno níže.

(Komponenty kromě Internal Stain, External Stain a FC Paste Stain mohou být smíchány s destilovanou vodou.)

Komponenta	Správný liquid k použití				
	CERABIEN™ZR FORMING LIQUID 100mL	Magic Formaer 20mL	MEISTER LIQUID 100mL	IS LIQUID 10mL/100mL	ES LIQUID 10mL/100mL
Margin	●	●			
Shade Base	●		●		
Opacious Body	●		●		
Cervical	●		●		
Body	●		●		
Enamel	●		●		
Translucent	●		●		
Luster	●		●		
Clear Cervical	●		●		
Mamelon	●		●		
Tissue	●		●		
Modifier	●		●		
Add-on	●		●		
Internal Stain				●	
External Stain					●
FC Paste Stain					●